

LAMINO HEPA-IV

Amino Acids Infusion 8% w/v; 200 ml and 500 ml Glass Bottle

LAMINO HEPA-IV

Amino Acids Infusion 8% w/v; 200 ml and 500 ml Glass Bottle

DESCRIPTION:

LAMINO HEPA-IV is a clear, colorless infusion, containing well balanced mixture of pure crystalline, essential, semi-essential and nonessential amino acids, which are in optimum concentration for protein synthesis.

INDICATIONS:

- Indicated as an amino acid (nitrogen) source in parenteral nutrition regimens. This use is appropriate when the enteral route is inadvisable, inadequate or not possible, as when:
 - Gastrointestinal absorption is impaired by obstruction, inflammatory disease or its complications, or antineoplastic therapy;
 - 🖔 Bowel rest is needed because of gastrointestinal surgery or its complications such as ileus, fistula or anastomotic leaks;
 - 🖔 Tube feeding methods alone cannot provide adequate nutrition.
- For treatment of hepatic encephalopathy in patients with cirrhosis/hepatitis & hepatic coma.

BENEFITS & FEATURES:

- Rich in branched chain amino acids (BCAA)
- Adequate protein supply during hepatic insufficiency
- Administration through peripheral vein

DOSAGE:

- Recommended daily dosage: For intravenous infusion 1.2 to 1.6 ml/kg body weight/hour.
- Dosage should be adjusted according to amino acid requirements depending upon age, body weight and condition of
 patients. Adequate calories should be administered simultaneously either orally or parenterally.

INFUSION RATE:

LAMINO HEPA-IV should be infused at 28 to 38 drops / minute for 70 kg body weight patient.

ADMINISTRATION:

LAMINO HEPA-IV is available as a sterile, non-pyrogenic, single dose container that can be administered through peripheral veins or by central venous route using non-pyrogenic I.V. administration set using aseptic technique.

PRESENTATION:

LAMINO HEPA-IV is available in 200 ml and 500 ml glass bottle.

La Renon Healthcare Private Limited

207-208 Iscon Elegance, Circle P, Prahlad Nagar Cross Roads, S.G. Highway, Ahmedabad-380015, Gujarat, India.

Phone: + 91-79-6616-8998, 2693-6656

E-mail: info@larenon.com | Web: www.larenon.com

